



# Strategies to Accelerate Process Development from Preclinical to Manufacturing for Gene Therapy

Dogan Ornek, Xingcheng Chen, Hassan Sakhtah, Chih-Cheng Chen, Nida S. Zubairy, Tao Xiang, Jared Kaufman, Linda Ma, Acelya Babali, Chi-Ming Yu, W. Roy Lin

Boston Institute of Biotechnology, LLC, 225 Turnpike Road, Southborough, MA 01772

Contact Information: [don.ornek@bostonbib.com](mailto:don.ornek@bostonbib.com) | +1-833-818-CDMO (2366) | [www.bostonbib.com](http://www.bostonbib.com)

## ABSTRACT

Gene therapy is a fast-growing industry. Demand on its products such as pDNA, AAV, retrovirus, lentivirus, etc. has significantly increased during last 10 years while its production challenges remain unresolved. It is a complex process with a number of unoptimized and low yielding steps. A typical process development (PD) cycle from preclinical to late phase manufacturing (Mfg) takes 3 to 5 years. It is desired to rapidly advance to manufacturing for early clinical candidates using well established and high yielding processes. The following strategies were successfully employed by BIB to accelerate the process development activities:

- (1) A platform approach for early stage process development,
- (2) Efficient DoE application for process optimization and characterization,
- (3) Interaction improvement between PD and Mfg to streamline scale-up, tech transfer and manufacture support,
- (4) Development of a new PAT tools application to understand process design space and hence develop a control strategy.

Case studies for BIB AAV and pDNA production platforms will be covered in this presentation.

## GENE THERAPY IS A COMPLEX PROCESS: PROCESS FLOW CHART FOR AAV PRODUCTION

| Plasmid Manufacturing           | Viral Manufacturing                                     |  | Primary Recovery                             | Downstream Processing                                | Analytical Support                             |
|---------------------------------|---|--|--|--|--|
| Vial thaw (E.coli)              | Vial thaw   |  | Harvest transfected cells                    | 2-step purification                                  | Identity: SDS-PAGE, MS, PCR                    |
| Fermentation                    | Adherent  | Suspension   | Cell lysis (mech., chem.)                    | Buffer exchange by TFF                               | Titer: qPCR, ddPCR, A260/A280                  |
| Plasmid recovery                | Cell Expansion 96/flat, roller bottle, multi-layer tray | Cell Expansion 5L wave or stirred-tank bioreactors | Nuclease treatment                           | Formulation and characterization                     | Potency: PFU, FFU                              |
| 2-step purification             | Serial passaging of cells                               | Serial passaging of cells                          | Concentration and dilution of harvest by TFF |  | Purity: TEM, MS, HCP-DNA                       |
| Formulation                     | Production/ Multi-layer tray, packed bed                | Production/ Stirred-tank, wave, perfusion          |  |  | Safety: Sterility, Endotox, Myco.              |
|                                 |   |  |  |  | Stability: pH, Osmo, Aggregate                 |
| Analytical Support              |   |  |  |  |  |
| Identity: Capsid ELISA, MS, PCR | Titer: qPCR, ddPCR                                      | Potency: TCID50, Gene Function Test                | Purity: SDS-PAGE, A260/A280, TEM, HCP-DNA    | Safety: Sterility, Endotox, Myco, RCV, Adenov. Agent | Stability: pH, Osmo, Aggregate, titer, potency |

## STRATEGIES TO ACCELERATE PROCESS DEVELOPMENT FROM PRECLINICAL TO MANUFACTURING


1. Development of platform for early-stage Process Development (PD)
2. Use of efficient DoE for late-stage PD
3. Improvement of interactions between PD and Mfg
4. Application of new Process Analytical Technology (PAT)

Note: Only BIB's AAV and pDNA production platforms are covered in this presentation.

## DEVELOPMENT OF PLATFORM TO ACCELERATE PROCESS DEVELOPMENT


### Strategies

1. Advanced expression toolbox with high-throughput screening
2. Standard media & feed solutions
3. Standard seed expansion & production conditions
4. Single use bioreactors and high-throughput analytics



## BIB's AAV MEDIA TOOLBOX PLATFORM FOR QUICK AND EFFICIENT PRODUCTION

BIB's platform is a collection of high performing media/reagents/cell lines for screening high AAV titer producing combos




- Identify top media/reagent/Cell combos based on Viral Genome (VG) production, growth, metabolism, and product quality

- Optimize AAV production in top performing combos using DoE
- Vector ratio, reagent: vector ratio, cocktail volume, infection VGD, Infection Duration

- Bioreactor test of top 1-3 media/cell combos
- Choose final media/cell combo based on growth, metabolism, GC and product quality

## BIB's pDNA PLATFORM FOR QUICK AND EFFICIENT PRODUCTION

BIB's pDNA platform includes strain development, fermentation, recovery, purification, and analytical methods. It is scalable to a thousand liter.




- BIB's platform is compatible with common industry strains such as DH5α & DH1

- Unique fed-batch high-cell-density fermentation reliably produces ≥200 mg/L pDNA (for typical ~6kb pDNA vaccines).

- Novel efficient primary recovery step that improves yield
- A single novel filtration step replaces the conventional 3-step primary recovery (flocculation → centrifugation → filtration) which reduces processing time & improves yield


- High-yielding purification process for isolating high purity super-coiled pDNA
- Overall downstream yields is >50% (overall yield 100 mg/L pDNA, >98% purity)

## CASE STUDY 1: EFFECTS OF MEDIA TOOLBOX ON AAV PRODUCTION




BIB platform resulted in 50-fold titer increase compared to the client's initial process

## CASE STUDY 2. BIB'S pDNA FERMENTATION IS SIMPLE AND SCALABLE



## CASE STUDY 2. BIB'S pDNA FERMENTATION PROVIDES CONSISTENT pDNA PRODUCTION



## SUMMARY FOR BIB PLATFORM APPROACH FOR GENE THERAPY

BIB successfully employed platform approach to accelerate gene therapy process development. It provides

- Timeline reduction
- Cost saving
- Robust process
- High yielding titer
- Reproducible product titer and quality